



SINGLE USE BIOPSY FORCEPS

Manufactured with the best Quality Medical Grade Materials



High quality flexible body for smooth insertion thus minimising the risk of scope damage

Ergonomic Handle designed for safe & easy operations, which precisely relay tiny finger movements to the distal jaws.



MIM serrated smooth Jaws. Available with & without Spike.
Wide opened jaws for sampling of large & diagnostically conclusive Biopsies.

VARIANTS:

| CODE. NO. | SHEATH/JAW O.D. | LENGTH | CUP TYPE |
|-----------|-------------------------|--------|---------------|
| ABF7-GS | 2.3mm/7Fr | 170cm | With Spike |
| ABF7-CS | 2.3mm/7Fr | 230cm | With Spike |
| ABF7-BS | 2.3mm/7Fr | 100cm | With Spike |
| ABF7-G | 2.3mm/7Fr | 170cm | Without Spike |
| ABF7-C | 2.3mm/7Fr | 230cm | Without Spike |
| ABF7-B | 2.3mm/7Fr | 100cm | Without Spike |
| ABF7-G-F | 2.3 mm/ 7Fr Fenestrated | 170 cm | Without Spike |
| ABF7-C-F | 2.3 mm/ 7Fr Fenestrated | 230 cm | Without Spike |
| ABF-B-F | 2.3 mm/ 7Fr Fenestrated | 100 cm | Without Spike |

DEVICE DESCRIPTION :

Glowtronics single use Biopsy Forceps are designed for endoscopic use to retrieve tissue samples for analysis. Forceps outer diameters are 2.3 mm and working lengths vary from 100cm to 230cm with serrated smooth jaws with and without Spike.

INTENDED USE :

This device is used to collect tissue or foreign body retrieval, from gastrointestinal tract using a Flexible endoscope.

WARNING :

Do not use this device for any purpose other than stated intended use.
If package is opened or damaged when received, do not use. Do not use if labeling is incomplete or illegible, In case Biopsy Forceps is supplied Unsterile – Ensure sterilization of the device before use

CAUTION :

This device is designed for single use only.
Read the entire IFU before use.
Refer to package label for minimum channel size required for this device.



Glowtronics Pvt Ltd
1D, Hootagalli Industrial Area,
Mysore – 570 018, India Tel: +91-821-240-2461/2462
Email: marketing@glowtronics.com Website: www.glowtronics.com

