

# SINGLE USE LAPAROSCOPIC SCISSORS

Manufactured with the best Quality Medical Grade Materials



**Hardened Precision ground / Corrosion Resistance Stainless Steel Blades for optimum performance.**  
**Dual- Edge smooth cutting action.**

**Wide opening to take in larger tissue for cutting, Allows Dissection and transection of fine tissue**  
**Does not loose its sharpness over extensive use**

**VARIANTS :**

CODE NO	DIAMETER	WORKING LENGTH	BLADE LENGTH
SLS3-18-C1	5mm	330mm	18mm
SLS4-18-C1	5mm	430mm	18mm

**DEVICE DESCRIPTION :**

Glowtronics disposable laparoscopic Scissors are intended for use in minimally invasive laparoscopic surgical procedures for Cut, dissect , specimen retrieval, coagulating, cauterizing and other manipulation of tissues and vessels during laparoscopic procedures. The insulated shaft is of 5mm diameter that connects to a handle with a standard monopolar male banana plug to be utilized for monopolar procedures when attached to their standard monopolar cables and their generators.

**INTENDED USE :**

This device is used to cut, dissect , coagulate tissue , retrieve foreign bodies, and remove tissues from within the gastrointestinal tract

**WARNING :**

This device is intended for Single - use only  
 Follow sterilization and other operating room standard procedures while using this device

**CAUTION :**

To avoid damage to patient, operator & instrument, become familiar with a specific instrument and its clamping or cutting mechanism prior to employing it in a surgical procedure.  
 Instrument is designed for grasping and holding soft tissue.  
 Attempting to cut staples or clips will damage the instrument.



**Glowtronics Pvt Ltd**  
 1D, Hootagalli Industrial Area,  
 Mysore – 570 018, India Tel: +91-821-240-2461/2462  
 Email: [marketing@glowtronics.com](mailto:marketing@glowtronics.com) Website: [www.glowtronics.com](http://www.glowtronics.com)