


# SINGLE USE LAPAROSCOPIC GRASPERS

Manufactured with  
the best Quality  
Medical Grade  
Materials



## VARIANTS:

CODE NO	DIAMETER	WORKING LENGTH	JAW LENGTH
SLG3-17-S1	5mm	330mm	17mm
SLG3-18-S1	5mm	330mm	18mm
SLG4-18-S1	5mm	430mm	18mm

## DEVICE DESCRIPTION :

Glowtronics Laparoscopic Graspers used to grasp tissue , retrieve foreign bodies, and remove tissues from within the gastrointestinal tract during endoscopic/laparoscopic surgical procedures. The insulated shaft is of 5mm diameter that connects to a handle with a standard monopolar male banana plug to be utilized for monopolar procedures when attached to their standard monopolar cables and their generators.

## INTENDED USE :

This device is used to grasp tissue , retrieve foreign bodies, and remove tissues from within the gastrointestinal tract

## WARNING :

This device is intended for Single - use only  
Follow sterilization and other operating room standard procedures while using this device

## CAUTION :

To avoid damage to patient, operator & instrument, become familiar with a specific instrument and its clamping or grasping mechanism prior to employing it in a surgical procedure.  
Instrument is designed for grasping and holding soft tissue.  
Attempting to cut staples or clips will damage the instrument.



Glowtronics Pvt Ltd  
1D, Hootagalli Industrial Area,

Mysore – 570 018, India Tel: +91-821-240-2461/2462

Email: [marketing@glowtronics.com](mailto:marketing@glowtronics.com) Website: [www.glowtronics.com](http://www.glowtronics.com)

