


## SINGLE USE BLADDER EVACUATOR STERILE / UNSTERILE



### VARIANTS:

BLE-01 (Soft Tubing Length : 40mm)

### DEVICE DESCRIPTION :

The Bladder Evacuator is made of semi-transparent materials. It is supplied in a sterile, peel-open package and is intended for one-time use only.

### INTENDED USE :

Bladder Evacuator is intended to be used for evacuating irrigation, tissue collection and/or straining functions during transurethral prostate or bladder surgery.

### WARNING :

All fittings to be tightened properly.

Select the correct nozzle end which is sized for the ACMI, STORZ, Olympus or Wolf Resectoscope.

The evacuator has to be assembled prior to use

Incase supplied Unsterile – Ensure sterilization of the device before use

If package is opened or damaged when received, do not use. Do not use if labeling is incomplete or illegible

### CAUTION :

These devices are for Single Use ONLY.

Check expiry date prior to use. Do not use after expiry date.

Do not Reuse, Re - process or Re-sterilize

This device must be stored in dry & cool conditions, away from extreme temperature.

Use of this device for any other purpose other than the stated intended use is neither advised nor recommended.



**Glowtronics Pvt Ltd**  
1D, Hootagalli Industrial Area,  
Mysore – 570 018, India Tel: +91-821-240-2461/2462  
Email: [marketing@glowtronics.com](mailto:marketing@glowtronics.com) Website: [www.glowtronics.com](http://www.glowtronics.com)



intertek



intertek

